



Certex[®]

Spinal Fixation System




Comprehensive Posterior Cervical System

Certex[®] Spinal Fixation System



The Certex OCT spinal fusion system offers a comprehensive solution for occipito-cervico-thoracic stabilization. The system features an innovative polyaxial screw design enabling high angulation in four positions.

- 4x40 degree favored angle provides 80 degrees of polyaxial motion with 62 degrees of standard motion
- Conical shape cap provides multiple point fixation to the rod while a recessed marker gives visual identification of the lead-in thread to minimize cross threading.
- Extensive implant selection and streamlined instrumentation provides the ability to address various surgical needs.



888.886.9354
cs@xtantmedical.com
xtantmedical.com
 664 Cruiser Ln,
 Belgrade, MT 59714

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.
 LIMITED WARRANTY and DISCLAIMER: Xtiant Medical products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed. WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.