


IRIX[®]-A


Lumbar Integrated Fusion System



Stand-alone device with
integrated locking mechanism

Irix[®]-A

Lumbar Integrated Fusion System




The Irix-A is an integrated stand-alone device for use in anterior lumbar interbody fusions. The zero profile with integrated zero-step locking mechanism provides 3 indicators of screw implantation.

➤ Integrated titanium ring encapsulated in PEEK provides increased torsional strength compared to a leading competitor.


➤ Self-tapping and self drilling screws are available in multiple lengths with 6 degrees of conical angulation.

➤ Frictional titanium plasma coating on bone interfacing surfaces aids in migration resistance, provides endplate visualization and aids in osteointegration.




Xray image of Irix-A showing zero profile.


Irix-A
33mm x 26mm Implants



Irix-A
38mm x 28mm Implants



Irix-A
42mm x 30mm Implants



☎ 888.886.9354
✉ cs@xtantmedical.com
🌐 xtantmedical.com

📍 664 Cruiser Ln,
Belgrade, MT 59714

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.

LIMITED WARRANTY and DISCLAIMER: Xtant Medical products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed. WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

Irix-A is a registered trademark of Xspine Systems, a subsidiary of Xtant Medical.
© Xtant Medical. All Rights Reserved. FM-C-MRK-31 (C) 05/2020