

SimpliMix™

Graft Delivery Device




Targeted Delivery Device

SimpliMix™

Graft Delivery Device

The SimpliMix™ Graft Delivery Device is designed for targeted delivery of hydrated allograft, autograft, or synthetic bone graft materials to an orthopedic surgical site, while maximizing material utilization.

- **Integrated Stylet** - Maximizes material utilization and simplifies delivery with an all-in-one configuration
- **Vented Configuration** - Vents air for better saturation and hydration of materials
- **Tapered Cap and Wide Bore Cannula** - Allows for the extrusion of high viscosity bone materials
- **Easy to use** - Available female luer adapter connects directly with a standard syringe for easy hydration, and large finger flanges allow for a better grip and is designed for one-handed ease-of-use



Item #	Description
400015	SimpliMix Graft Delivery Device 15cm



☎ 888.886.9354
 ✉ cs@xtantmedical.com
 🌐 xtantmedical.com
 📍 664 Cruiser Ln,
 Belgrade, MT 59714

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.
 LIMITED WARRANTY and DISCLAIMER: Xtant Medical products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed. WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

SimpliMix is a trademark of Xtant Medical.
 © Xtant Medical. All Rights Reserved. FM-C-MRK-124(A) 05/2021